

# Experimental and model-based investigation of twin screw granulation

*Ashish Kumar*

*Lab Meeting, FFW      19 October, 2015*

LABORATORY OF PHARMACEUTICAL PROCESS ANALYTICAL TECHNOLOGY


FACULTY OF PHARMACEUTICAL SCIENCES

BIOMATH, DEPARTMENT OF MATHEMATICAL MODELLING, STATISTICS AND BIOINFORMATICS


FACULTY OF BIOSCIENCE ENGINEERING

# Continuous manufacturing line

## Consigma™-25 system




# At appropriate time-scales and conditions, granulation is in steady state





## Two key implications

1. Fluxes are roughly constant (Dynamics are transient)
2. If feed is constant, product quality is consistent!



# Twin-Screw Granulator applies High Shear Wet Granulation



# Twin-screw granulation process development

Areas under study:

- Granulation time and mixing.
- Aggregation and breakage rates.

# Twin-screw granulation process development


Areas under study:

- Granulation time and mixing.
- Aggregation and breakage rates.

# Consigma™-1 system

(GEA pharma systems, Collette)

## Open barrel of a twin screw granulator







# Consigma™- 1 experiments

Lactose/PVP (97.5/2.5) premix was granulated with distilled water

Factors:


Parameters	Low	High
Throughput	10 Kg/h	25 Kg/h
Liquid-solid ratio	4.58 %	6.52%
Screw speed	500 RPM	900 RPM




Responses:

Particle characterization by Dynamic Image Analysis  
(Location 1, 3, 5)


# Average Feret diameter vs Aspect ratio



# Comparing average Feret diameter



# Granulation is result of particle population dynamics



# Population balance equation

$$\frac{\partial n(t, x)}{\partial t} = \frac{Q_{in}}{\tilde{V}} n_{in}(x) - \frac{Q_{out}}{\tilde{V}} n_{out}(x) \quad \text{GSD balance}$$

Aggregation  
term


$$+ \frac{1}{2} \int_0^x \beta(t, x - \varepsilon, \varepsilon) n(t, x - \varepsilon) n(t, \varepsilon) d\varepsilon \\ - n(t, x) \int_0^\infty \beta(t, x, \varepsilon) n(t, \varepsilon) d\varepsilon$$

Breakage  
term


$$+ \int_0^\infty b(x, \varepsilon) S(\varepsilon) n(t, \varepsilon) d\varepsilon \\ - S(x) n(t, x)$$

$\beta$  = aggregation rate  
 $S$  = selection rate  
 $b$  = breakage function


# Experimental and simulated data have a good agreement



# Particle population dynamics during granulation





# Experimental and simulated data have a good agreement






# Including effect of granulator design on granule size distribution




**IS IT EYECON RESULTS USEABLE?**

# Comparison of PSD analysis methods

- Sieve analysis
- 3D imaging techniques (Eyecon™)




# Three scenarios




# Process settings for the runs

	Run 1	Run 2	Run 3	Run 4
Throughput (kg/h)	25	10	10	25
Liquid addition (%)	9	9	8	9
Screw	1x6	1x6	2x6	2x6
Screw Speed (RPM)	500	900	900	900
Stagger angle (°)	60	60	60	60
Temperature (°C)	25	25	25	25
Binder Addition	Wet	Wet	Wet	Wet

# Difference between Eyecon before dryer and sieve after dryer




# Difference between Eyecon before dryer and sieve after dryer






# Difference between Eyecon before dryer and sieve before dryer






# Difference between Eyecon before dryer and sieve before dryer


More granules in eyecon measurement




More granules in sieve before dryer measurement



# Difference between sieve before dryer and sieve after dryer



# Difference between sieve before dryer and sieve after dryer



# Conclusions

Along with experimental study, **an improved insight** can be obtained by model-based analysis.

**High throughput processing can be achieved** by increasing the liquid-solid ratio and screw speed.

**PBM requires further development** to include screw geometry effect and validation.

Eyecon is **showing some promising results**, and can be used for studies in future.

# Aknowledgements

Prof. Thomas De Beer


Prof. Ingmar Nopens

Prof. Krist V. Gernaey

Jurgen Vercruysse

Fien De Leersnyder





Laboratory of Pharmaceutical Process Analytical Technology

Ashish.Kumar@UGent.be

- Extra slides