

# Linking granulation performance with residence time & liquid distributions in twin-screw granulation


*Ashish Kumar*

*LPPAT-lab meeting*




# Traditional to new granulation method

High-shear mixer




Batch

Drum




From minutes to hours

Fluidised-bed





Continuous

Twin-screw




# At appropriate conditions, granulation is in steady state



**Steady state**

*transfer in*  $\approx$  constant  $\approx$  *transfer out*


$$\frac{d[P_m]}{dt} \approx 0 \approx \frac{d[G_m]}{dt}$$



## Two key implications



1. Fluxes are roughly constant (Dynamics are transient)
2. Same amount of time is to complete all sub-processes

# High Shear Wet Granulation involves different rate processes




# Having many time-scales is challenging

## Residence time-scale



# Having many time-scales is challenging

## Mixing time scale




# Residence time and moisture distributions effect on the granulation performance

## Screw Configuration

- Number of kneading discs (4, 6, 2x6)
- Stagger angle (30°, 60°, **120°**)

## Process parameters

- Material throughput (10-25 kg/h)
- Screw speed (500-900 rpm)
- Liquid-to-solid ratio (6-8%)



spike

Fines  
< 150  $\mu\text{m}$

Yield fraction  
> 150 to <1400  $\mu\text{m}$

Oversized  
> 1400  $\mu\text{m}$

# Analysis of distributions in twin-screw granulation


Measurement by distributions

Results


Summary




# Tracer concentration in granules measured by NIR chemical imaging



# Tracer maps used to measure distributions



# Qualitative assessment of the RTD profiles



$$\tau = \frac{\int_0^{\infty} t \cdot e(t) dt}{\int_0^{\infty} e(t) dt}$$

Mean residence time ,  $\tau$   
(a measure of the mean of the distribution)


$$\sigma^2 = \frac{\int_0^{\infty} (t-\tau)^2 \cdot e(t) dt}{\int_0^{\infty} e(t) dt}$$

Variance,  $\sigma^2$   
(width of the distribution)

$$Pe = \frac{UL}{D}$$

Péclet Number,  $Pe$   
 $\left( \frac{\text{Rate of axial transport by convection}}{\text{Rate of axial transport by dispersion}} \right)$

# Qualitative assessment of the moisture maps




Shannon Entropy based Mixing Index

$$H(X) = \sum_{j=1}^n P(X_j) \log_{200} (1 / P(X_j))$$

$$MI = -\frac{1}{\log_{200}(n)} \sum_{j=1}^n P(X_j) \log_{200} P(X_j)$$

FFT to obtain Frequency and amplitude



# Analysis of distributions in twin-screw granulation


RTD Measurement by Chemical Imaging

Results


Summary

# Increase in L/S lubricates moving parts but flow is sluggish

## Torque




## Mean residence time




Torque (N=42; DF=32; R2=0.95), Mean residence time (N=41; DF=30; R2=0.97) Conf. lev.=0.95

# Mixed-flow transport at a high screw speed & plug-flow transport when more kneading discs

## Normalised variance



## Péclet number



Normalised variance (N=46; DF=39; R2=0.87), Péclet number (N=43; DF=35; R2=0.84) Conf. lev.=0.95


Axial mixing  $\uparrow$  when Normalised variance  $\uparrow$  and Péclet number  $\downarrow$




# Throughput force is more dominant at a high screw speed

● MFR (low)  
■ MFR (high)

**Mean residence time**




**Peclet number**






# Increase no. of kneading and in L/S led to improved liquid distribution in bulk






Mixing index



Frequency


Amplitude



Mixing index (N=43; DF=34; R<sup>2</sup>=0.80), Frequency (N=41; DF=32; R<sup>2</sup>=0.84), Amplitude (N=45; DF=38; R<sup>2</sup>=0.88) Conf. lev.=0.95


S/L mixing  $\uparrow$  when MI  $\uparrow$ , Frequency  $\downarrow$  and Amplitude  $\downarrow$

# Trade-off existed between high throughput and number of kneading discs




# Increase in L/S and no. of kneading led to improved liquid distribution, hence less fines

## Fine




## Yield

## Oversized



Fine (N=49; DF=39; R2=0.93), Yield (N=38; DF=28; R2=0.84), Oversized (N=47; DF=38; R2=0.92) Conf. lev.=0.95



Fine (N=49; DF=39; R2=0.93), Yield (N=38; DF=28; R2=0.84), Oversized (N=47; DF=38; R2=0.92) Conf. lev.=0.95

# Analysis of distributions in twin-screw granulation

RTD Measurement by Chemical Imaging

Results

Summary

# The results showed that..

..material throughput and number of kneading discs dictate solid-liquid mixing.

...till good mixing kneading discs are not there TSG should better be operated at lower throughput.

.. non-conventional screw elements with modified geometries should be explored for improvement in solid-liquid mixing.

Prof. Thomas De Beer  
Prof. Ingmar Nopens  
Prof. Krist V. Gernaey

Valérie Vanhoorne  
Fien De Leersnyder  
Dr. Jurgen Vercruysse

Maunu Toiviainen  
Dr. Mikko Juuti

Maija Alakarjula  
Prof. Jarkko Ketolainen



UNIVERSITY OF  
EASTERN FINLAND

Discrete element method (DEM):  
detailed simulation of wet granulation



# Involved physics in granulation

- Mixing
- Wetting
- Aggregation (PBM)
- Breakage (PBM)


# Challenges of PBM

- Does not include geometrical effect
- Rate of sub-processes (i.e. aggregation, breakage) is lumped into kernels which are
  - Size dependent
  - Energy/shear dependent
- Need extensive experimental calibration


# Every contact and neighborhood matter!



# DEM for detailed investigation




# Particle scale DEM





# Initial parameter values

Parameter	values
Throughput	10-25 kg/h
Material density	1500 kg/m <sup>3</sup>
Bulk density	500-750 kg/m <sup>3</sup>
Particle diameter	30-70 $\mu\text{m}$
Number of particles	~60- 150 million # /sec

# A lot of detail can be extracted





# Sheared-box simulation





# Breakage is dominant mechanism at high fill ratio





Laboratory of Pharmaceutical Process Analytical Technology

BIO MATH

Model-based analysis and optimization of bioprocesses

Ashish.Kumar@UGent.be