

Experimental and model-based investigation of twin screw granulation

Ashish Kumar

Lab Meeting, FFW 18 August 2014

LABORATORY OF PHARMACEUTICAL PROCESS ANALYTICAL TECHNOLOGY


FACULTY OF PHARMACEUTICAL SCIENCES

BIOMATH, DEPARTMENT OF MATHEMATICAL MODELLING, STATISTICS AND BIOINFORMATICS


FACULTY OF BIOSCIENCE ENGINEERING

Continuous manufacturing line

Consigma™-25 system




At appropriate time-scales and conditions, granulation is in steady state





Two key implications

1. Fluxes are roughly constant (Dynamics are transient)
2. If feed is constant, product quality is consistent!



Twin-Screw Granulator applies High Shear Wet Granulation



Twin-screw granulation process development

Areas under study:

- Granulation time and mixing.
- Aggregation and breakage rates.

Twin-screw granulation process development


Areas under study:

- Granulation time and mixing.
- Aggregation and breakage rates.

Consigma™-1 system

(GEA pharma systems, Collette)

Open barrel of a twin screw granulator




Consigma™- 1 experiments

Lactose/PVP (97.5/2.5) premix was granulated with distilled water

Factors:

Parameters	Low	High
Throughput	10 Kg/h	25 Kg/h
Liquid-solid ratio	4.58 %	6.52%
Screw speed	500 RPM	900 RPM



1

2

kneading block 1

3

4


kneading block 2

5


Responses:

Particle characterization by Dynamic Image Analysis
(Location 1, 3, 5)


Average Feret diameter vs Aspect ratio



Comparing average Feret diameter



Granulation is result of particle population dynamics



Population balance equation

$$\frac{\partial n(t, x)}{\partial t} = \frac{Q_{in}}{\tilde{V}} n_{in}(x) - \frac{Q_{out}}{\tilde{V}} n_{out}(x) \quad \text{GSD balance}$$

Aggregation
term


$$+ \frac{1}{2} \int_0^x \beta(t, x - \varepsilon, \varepsilon) n(t, x - \varepsilon) n(t, \varepsilon) d\varepsilon \\ - n(t, x) \int_0^\infty \beta(t, x, \varepsilon) n(t, \varepsilon) d\varepsilon$$

Breakage
term


$$+ \int_0^\infty b(x, \varepsilon) S(\varepsilon) n(t, \varepsilon) d\varepsilon \\ - S(x) n(t, x)$$

β = aggregation rate
 S = selection rate
 b = breakage function


Experimental and simulated data have a good agreement





Particle population dynamics during granulation



Experimental and simulated data have a good agreement




Including effect of granulator design on granule size distribution




IS IT EYECON RESULTS USEABLE?

Comparison of PSD analysis methods

- Sieve analysis
- 3D imaging techniques (Eyecon™)




Three scenarios




Process settings for the runs

	Run 1	Run 2	Run 3	Run 4
Throughput (kg/h)	25	10	10	25
Liquid addition (%)	9	9	8	9
Screw	1x6	1x6	2x6	2x6
Screw Speed (RPM)	500	900	900	900
Stagger angle (°)	60	60	60	60
Temperature (°C)	25	25	25	25
Binder Addition	Wet	Wet	Wet	Wet


Difference between Eyecon before dryer and sieve after dryer



Difference between Eyecon before dryer and sieve after dryer




Difference between Eyecon before dryer and sieve before dryer




Difference between Eyecon before dryer and sieve before dryer


More granules in eyecon measurement




More granules in sieve before dryer measurement



Difference between sieve before dryer and sieve after dryer



Difference between sieve before dryer and sieve after dryer



Conclusions

Along with experimental study, **an improved insight** can be obtained by model-based analysis.

High throughput processing can be achieved by increasing the liquid-solid ratio and screw speed.

PBM requires further development to include screw geometry effect and validation.

Eyecon is **showing some promising results**, and can be used for studies in future.

Aknowledgements

Prof. Thomas De Beer


Prof. Ingmar Nopens

Prof. Krist V. Gernaey

Jurgen Vercruysse

Fien De Leersnyder





Laboratory of Pharmaceutical Process Analytical Technology

Ashish.Kumar@UGent.be