
GEA Event

Ashish Kumar

At appropriate time-scales and conditions, granulation is in steady state




 $transfer in \approx constant \approx transfer out$

$$\frac{d[P_m]}{dt} \approx 0 \approx \frac{d[G_m]}{dt}$$


Two key implications

- 1. Fluxes are roughly constant (Dynamics are transient)
- 2. If feed is constant, product quality is consistent!

Scope: Desired Outcome



Twin-Screw Granulator applies High Shear Wet Granulation




Granule Size Distribution

Key questions for twin-screw granulation process development

What affects granulation time and mixing?

What control degree of rate processes involved in desired quality of granules?

Residence time distribution to know the granulation time and mixing



Residence time distribution to know the granulation time and mixing


Screw Configuration

- Number of kneading discs
- Stagger angle




Process parameters


- Material throughput
- Screw speed




Tracer concentration in granules produced was measured using NIR chemical imaging



API Map was used to measure RTD



Conceptual model to include three main components of RTD




Modified Tank-in-Series model used

Modified Tank-In-Series model

$$e(\theta) = \frac{b[b(\theta - p)]^{n-1}}{(n-1)!} exp^{[-b(\theta - p)]}$$

where,
$$b = \frac{n}{(1-p)(1-d)}$$



Analysis of residence time distribution in twin-screw granulation


RTD measurement by Chemical imaging **Model Formulation**

Outcomes

Measurement based - Mean residence time

Number of CSTR - Plug flow fraction
- Dead volume fraction

Dead volume fraction



$$\tau = \frac{\int_0^\infty t \cdot e(t)dt}{\int_0^\infty e(t)dt}$$

Mean residence time, τ (a measure of the mean of the distribution)

$$\sigma^2 = \frac{\int_0^\infty (t - \tau)^2 \cdot e(t) dt}{\int_0^\infty e(t) dt}$$

Variance, σ^2 (width of the distribution)

Analysis of residence time distribution in twin-screw granulation


Chemical imaging based RTD measurement Model Formulation

Outcomes


Measurement based - Mean residence time

Number of CSTR
 Plug flow fraction
 Dead volume fraction


Residence time reduces with increase in Measure of the mean of the distribution screw speed




Residence time reduces with increase in throughput...but not always




Residence time increases with increase in number of kneading discs.



Residence time reduces with increase in stagger angle.



Mean of the residence time distribution



Analysis of residence time distribution in twin-screw granulation

Chemical imaging based RTD measurement

Model Formulation


Outcomes

Measurement based - Mean residence time
- Mean centred variance


- Number of CSTR
- Plug flow fraction
- Dead volume fraction

Parameters of the TIS model estimated using experimental RTD based on least SSE


$$e(\theta) = \frac{b[b(\theta - p)]^{n-1}}{(n-1)!} e^{[-b(\theta - p)]}$$




Plug flow component of the RTD




Plug flow fraction decreases with increase in screw speed and throughput




Mixed flow component of the RTD




Material throughput controls mixing which reduces with increase in throughput



Mixed flow component of the RTD



Dead zone increases with screw speed, and reduces with increase in kneading discs



RTD analysis showed that

Screw speed controls the residence time,

Material throughput controls mixing.

Key questions for twin-screw granulation process development

What affects granulation time and mixing?

What control degree of rate processes involved in desired quality of granules?

Consigma™-1 system (GEA pharma systems, Collette)

Open barrel of a twin screw granulator




Consigma - 1 experiments

Lactose/PVP (97.5/2.5) premix was granulated with distilled water


Factors:

| Parameters | Low | High |
|--------------------|---------|---------|
| Throughput | 10 Kg/h | 25 Kg/h |
| Liquid-solid ratio | 4.58 % | 6.52% |
| Screw speed | 500 RPM | 900 RPM |




Particle characterization by Laser Obscuration Time technique Responses: (Location 1, 3, 5)

Comparing average Feret diameter




Comparing average Feret diameter

At high Throughput and L/S



Average Feret diameter vs Aspect ratio



Population balance equation

$$\left| \frac{\partial n}{\partial t} + \nabla \cdot \left(\mathbf{u}_{p} n \right) - \nabla \cdot \left(D_{p} \nabla n \right) \right| =$$

Granule size distribution

$$N \delta(V - v_o) - \frac{\partial}{\partial V} (G n)$$

Nucleation and growth term

$$+\frac{1}{2}\int_{0}^{V}\beta(v,V-v)n(v)n(V-v)dv-n(V)\int_{0}^{\infty}\beta(v,V)n(v)dv$$
Aggregation term

$$+\int_{V}^{\infty} \Gamma(\Phi)b(V;\Phi)n(\Phi)d\Phi - \Gamma(V)n(V)$$

Breakage term

N = nucleation rate

G = growth rate


 β = aggregation rate

 Γ = breakage rate


b = daughter distribution

 v_0 = nuclei size

Ignoring effect of granulator design on granule size distribution



Multiple factors contribute to the variation in granule size distribution



Conclusions

Along with experimental study, an improved insight can be obtained by model-based analysis (as we already seen in case of RTD).

Screw speed controls the residence time, while the material throughput controls mixing.

High throughput can easily be achieved by simultaneously increasing the feed rate and screw speed.

What next...

Together with a Granule size distribution study it will be confirmed which mixing regime is most desirable.

In further study we will investigate material properties influence on the RTD and mixing.

Utilise the mixing and residence time information for mechanistic modeling of the TSG.

Thomas De Beer Ingmar Nopens Krist V. Gernaey UNIVERSITEIT GENT

Jurgen Vercruysse Valérie Vanhoorne

Maunu Toiviainen Panouillot Pierre-Emmanuel Mikko Juuti

Kris Schoeters





Laboratory of Pharmaceutical Process Analytical Technology



Ashish.Kumar@UGent.be